

Soft contact lenses

appenzeller
kontaktlinsen®
passt.

i-SOFT

3-monthly soft contact lenses



My yearly supply of soft
contact lenses. Per eye.

info@appenzeller-kontaktlinsen.ch
www.appenzeller-kontaktlinsen.ch



i-SOFT

passt.

3-monthly soft contact lenses

See clearly and distinctly. Nothing more. But nothing less either. This is guaranteed with the i-SOFT by Appenzeller Kontaktlinsen. And that for three months. Not more. And not less either.

If the contact lens is lost or damaged, its replacement is less expensive than for a 6-monthly or a 12-monthly lens. Ideal for patients who don't want to take any risks and don't want to be without quality contact lenses under any circumstances.

→ i-SOFT Family

i-SOFT
i-SOFT TD / TP
i-SOFT proASSIST
i-SOFT proASSIST TD / TP
i-SOFT Multi
i-SOFT Multi TD / TP

→ Exchange frequency

3 months

→ Front surface design

Spherical	
Prism ballast	Marker in 0° and 180°
Dynamic stabilisation	Marker in 0° and 180°

→ Back surface design

- Spherical or toric
- Bi-curve, available in three different steps of peripheral flattening according to the corneo-scleral profile

Description	Design
1	low peripheral flattening
2	normal peripheral flattening (standard)
3	high peripheral flattening



→ Product range

Base curve (diameter depending)	7.5 mm to 9.2 mm in 0.1 mm steps
Diameter	13.3 mm / 13.8 mm / 14.3 mm / 14.8 mm / 15.3 mm
Flattening factor «AF»	1 / 2 / 3
Spherical power range	± 30.0 D in 0.25 D steps
Cylinder (for toric versions)	-0.5 D to -8.0 D in 0.25 D steps
Axis (for toric versions)	0° to 180° in 1° steps
Addition (for multifocal versions)	0.5 D to 3.5 D in 0.25 D steps
Risk-factor (for proASSIST versions)	Moderate / High / Very High

→ Materials

Material	Tint
Unisil Silicon Hydrogel 50	Clear / Blue
GM Advance 49%	Clear / Blue
GM Advance 58%	Clear / Blue
Contaflex 75%	Clear

→ Exchange/return periods

Date of the delivery note +1 month	Refund minus an excess
Date of the delivery note > 1 month to 1.5 months	50% Refund

→ Guarantee

Date of the delivery note +1 month	100% Refund on damage
------------------------------------	-----------------------

→ Fitting recommendations

Material																		
Unisil Silicon Hydrogel 50*																		
GM Advance 49%																		
GM Advance 58%																		
Contaflex 75%																		
Corneal radii mean value <small>(r1 + r2) / 2</small>		7.20	7.30	7.40	7.50	7.60	7.70	7.80	7.90	8.00	8.10	8.20	8.30	8.40	8.50	8.60	AF	
Cornea Ø (mm)	Small 11.0-11.4	7.70	7.80	7.90	8.00	8.10	x	x	x	x	x	x	x	x	x	x	AF 2	
	Medium 11.5-12.1	7.70	7.80	7.90	8.00	8.10	8.20	8.30	8.40	8.50	8.50	8.60	8.70	x	x	x	AF 2	
	Large >12.1	x	x	x	8.00	8.10	8.20	8.30	8.40	8.50	8.60	8.70	8.80	8.90	9.00	9.10	AF 2	

Diameter i-SOFT: Ø Cornea (HVID) +2.5 mm

*Unisil 50: Base curve 0.1 to 0.2 mm flatter than indicated in the table.



passt.

COMFORT AND SAFETY

Our recommendation for the care regime of individual soft contact lenses:

Individual products or sets.

Compiled to your wishes.

Contact lenses are not just a product. Contact lenses are a process that begins with your consultation. And continues after a successful fitting. Because contact lenses need to be inserted and removed. And because they need to be cleaned and stored. That's why we don't just supply you with contact lenses made to order, but also with know-how made to order. This includes our Professional Service, which is there for you in word and action. It also includes care products and accessories that give your patients the greatest possible feeling of comfort, health and safety.

Clean and rinse, disinfect and store : This is the sequence. And accordingly, we offer you products that are precisely matched to each other in exactly the combination you want. Put together nicely as a set wrapped in thin transparent film, which can be removed without scissors.

Combination example:

- 1x Cleaner plus 1x Saline plus 1x Peroxid plus 1x Enzyme
- 1x Cleaner plus 1x Saline plus 1x All-in-One Soft

See also brochure Appenzeller Kontaktlinsen Comfort and safety

1. Cleaning Daily

Appenzeller Kontaktlinsen® Cleaner

Alcohol-based cleaner

For daily care of rigid gas-permeable contact lenses (also GP-plasma-coated) and soft contact lenses.

2. Rinsing Daily

Appenzeller Kontaktlinsen® Saline

Saline as inserting solution

For all types of contact lenses - except scleral lenses - and for rinsing of soft and rigid gas-permeable contact lenses. Both after cleaning the lens and before putting the lens on the eye.



3. Disinfection and storage Daily

Appenzeller Kontaktlinsen® Peroxyd

Disinfection and neutralisation in 1 hour

The disinfection and neutralisation can be carried out very simply, quickly and a 100% efficiently. Appenzeller Kontaktlinsen make it possible with this peroxide system. For conventional soft lenses, disposable lenses and rigid gas-permeable lenses (RGP)

Appenzeller Kontaktlinsen® All-in-One Soft

Disinfectant

For all those who do not opt for the peroxide system and its unusually rapid effect (see column on the left): All-in-One Soft – the multipurpose solution with sodium hyaluronate for longer effect, longer wearing time, longer feeling of comfort.

4. Protein removal Once a week

Appenzeller Kontaktlinsen® Enzyme

Intensive cleaning

Protein deposits build up on the lens every day. In addition to the daily neutralisation of the lenses, a weekly protein removal is needed when using the peroxide system.





Little bowls of plastic...

passt.

It was Mr Feinbloom. The optician was the first one to use PMMA to manufacture contact lenses. That was in the 1930's. A big window to the world opened to millions of people since, thanks to those little bowls of plastic.

We put a great deal of effort into Research & Development to ensure that our contact lenses always meet the latest technological possibilities.

Our Professional Services team advises you competently, personally and patiently. They will answer every question and make sure you get the right contact lens in any case.

i-SOFT



3-monthly soft contact lenses
My yearly supply of soft contact lenses.
Per eye.

Personnelle



Unifocale – Multifocale – Style de vie numérique
My contact lens. As unique as my fingerprint.

Excellent AS Excellent MK



Single vision – Multifocal – Bifocal
Classic, practical, excellent:
You can count on this contact lens.

i-MAP AS



Single vision – Multifocal – Bifocal
High-tech for my everyday life:
Coping with any situation.

... big window to the world.

i-NIGHT



Orthokeratology at the highest level
Wear overnight – see during the day:
Without spectacles, without contact lenses.

Personnelle-KK Excellent-KK



Keratoconus contact lenses
High grade of keratoconus,
high demand: fits like a mould.

i-MAP



For fitting with irregular cornea
Fits when nothing else will.

i-MATRIX



Scleral lenses
Corneal damage, keratoconus,
keratoplasty? Still fits.

*pro*ASSIST



Professional support with progressive myopia
Smart contact lenses instead of thick
spectacle lenses.



passt.

Appenzeller Kontaktlinsen AG
Hauptstrasse 22
CH-9042 Speicher AR
Telephone +41 71 344 20 00

Appenzeller Kontaktlinsen AG
Kemptener Strasse 8
D-88131 Lindau
Telephone +49 69 332 962 31

Appenzeller Kontaktlinsen AG
Postfach 29
A-6973 Höchst
Telephone +43 72 088 16 97

info@appenzeller-kontaktlinsen.ch www.appenzeller-kontaktlinsen.ch

10 Reasons Why to choose Appenzeller Kontaktlinsen

→ Fitting success

We do everything we can to understand your wishes and ideas. So that your fitting is a complete success – for you and your patients.

→ Technological advantage

We put a great deal of effort into Research & Development to ensure that our contact lenses always meet the latest technological possibilities.

→ Swiss made

We demand the highest standards of the manufacturing quality of our contact lenses, so that they match your equally high expectations of Swiss made.

→ Reproducibility

We help you keeping your patients satisfied by ensuring that you receive the same lens as before with your repeat order.

→ Open ears

Our Professional Services team advises you competently, personally and patiently. They will answer every question and make sure you get the right contact lens in any case.

→ Custom made

Each contact lens is individually made to order for you, so that it matches your specifications a 100%.

→ Right to exchange

We send out our contact lenses with or without the right to exchange, so you can make the best choice for your needs. The option "with right to exchange" gives you financial security: you can order a different lens later. If you do so within the exchange period, you pay only a small excess. The date on the delivery note is decisive for the deadline.

- Soft contact lenses:
 - *i*-SOFT 3-monthly lenses: 1 month
 - 6-monthly lenses: 2 months
 - 12-monthly lenses: 3 months
- Rigid gas-permeable lenses: 3 months

→ Price advantage

We offer you favourable terms of payment. You can get extra benefits if you pay in advance. We will happily let you know about the current conditions on request.

→ Guarantee scheme

We take responsibility for our work so you can rely on us and our contact lenses with guarantee. Guarantee period:

- *i*-SOFT 3-monthly lenses: 1 month
- 6-monthly lenses: 2 months
- 12-monthly lenses: 3 months

Broken lens: Return us the lens within the guarantee period and we will replace it. We will credit it on your next invoice accordingly. The date on the delivery note is decisive for the guarantee period.

→ Long-term vision

We strive for a lasting business relationship with you, so that you can also count on us in the long term.

